
IPI549

产品编号: D60065

CAS: 1693758-51-8

分子式: C₃₀H₂₄N₈O₂

纯度: 99%

InChi Key: InChIKey=XUMALORDVCFWKV-IBGZPJMESA-N

Smiles: C[C@H](NC(=O)c1c(N)nn2cccnc12)c1cc2cccc(C#Cc3cnn(C)c3)c2c(=O)n1-c1cccc1

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PI3K γ ; PI3K

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Eganelisib (IPI-549) is an inhibitor of PI3K γ (IC₅₀ = 16, 3,200, 3,500, and >8,400 nM for PI3K γ , PI3K α , PI3K β , and PI3K δ , respectively).