
ARV-825

产品编号: D60075

CAS: 1818885-28-7

分子式: C₄₆H₄₇ClN₈O₉S

纯度: 97%

InChi Key: InChIKey=RWLOGRLTDKDANT-TYIYNFKSA-N

Smiles: CC=1C2=C(N3C([C@H](CC(NC4=CC=C(OCCOCCOCCOCCNC5=C6C(C(=O)N(C6=O)C7C(=O)NC(=O)CC7)=CC=C5)C=C4)=O)N=C2C8=CC=C(Cl)C=C8)=NN=C3C)SC1C

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: BRD4 BD2; BRD4 BD1; Epigenetic Reader Domain

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: ARV-825, a heterobifunctional PROTAC (Proteolysis Targeting Chimera) that recruits BRD4 to the E3 ubiquitin ligase cereblon leading to fast, efficient, and prolonged degradation of BRD4 in all BL cell lines tested.