
E3 Ligase Ligand-Linker Conjugates 2

产品编号: D60076

CAS: 1818885-63-0

分子式: C₂₇H₃₂N₄O₈

纯度: 98%

InChi Key: InChIKey=QTFDOBJYFQFDIO-UHFFFAOYSA-N

Smiles: Nc(cc1)ccc1OCCOCCOCCOCCNc1cccc(C(N2C(CCC(N3)=O)C3=O)=O)c1C2=O

外观: 固体粉末

作用靶点: Cereblon; Ligand for E3 Ligase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Pomalidomide-PEG4-Ph-NH₂ is a synthesized E3 ligase ligand-linker conjugate that incorporates the Pomalidomide based cereblon ligand and 4-unit PEG linker used in PROTAC technology.