

---

## Hemokinin 1(mouse)

产品编号: D60092

CAS: 208041-90-1

分子式:  $C_{61}H_{100}N_{22}O_{15}S$

纯度: 98%

外观: 固体粉末

作用靶点: Human NK1 receptor; Human NK2 receptor

溶解性: DMSO: 100 mg/mL (70.74 mM; Need ultrasonic)

保存条件: -20°C

产品介绍: Hemokinin 1 (mouse) is a selective agonist of neurokinin-1 receptor, with  $K_i$  of 0.175 nM and 560 nM for human NK1 receptor and human NK2 receptor, respectively.