

---

## DMA(trihydrochloride)

产品编号: D60093

CAS: 2095832-33-8

分子式:  $C_{27}H_{31}Cl_3N_6O_2$

纯度: 99%

外观: 固体粉末

作用靶点: HeLa cell; MCF7 cell; Others

溶解性: Soluble in H<sub>2</sub>O

保存条件: -20°C

产品介绍: DMA trihydrochloride is a fluorescent compound ( $\lambda_{ex}=340$  nm,  $\lambda_{em}=478$  nm).