
Tafluprost

产品编号: D60094

CAS: 209860-87-7

分子式: C₂₅H₃₄F₂O₅

纯度: 98%

InChi Key: InChIKey=WSNODXPBBALQOF-VEJSHDCNSA-N

Smiles: CC(C)OC(=O)CCC/C=C/[C@H]1[C@@H](O)C[C@@H](O)[C@@H]1\C=C\C(F)(F)COc1ccccc1

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: Prostaglandin Receptor

溶解性: DMSO:250 mg/mL (552.45 mM)

保存条件: 2-8°C

产品介绍: Tafluprost, a prostaglandin analog, is the selective agonist of fluoroprostaglandin (FP) receptor PGF₂α.