
DEGARELIX

产品编号: D60100

CAS: 214766-78-6

分子式: $C_{82}H_{103}ClN_{18}O_{16}$

纯度: 99%

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: GnRHR; GNRH Receptor

溶解性: DMSO: 10 mg/mL (6.13 mM; Need ultrasonic) H₂O: 5 mg/mL (3.06 mM; ultrasonic and warming and heat to 60°C)

保存条件: -20°C, 避光

产品介绍: Degarelix is a competitive and reversible gonadotropin-releasing hormone receptor (GnRHR/LHRHR) antagonist. Degarelix can be used for prostate cancer research