
MLCK INHIBITOR PEPTIDE 18

产品编号: D60107

CAS: 224579-74-2

分子式: C₆₀H₁₀₅N₂₃O₁₁

纯度: 98%

外观: 固体粉末

作用靶点: CaMK II

溶解性: H₂O: 100 mg/mL (75.49 mM; Need ultrasonic)

保存条件: -20°C

产品介绍: MLCK inhibitor peptide 18 is a myosin light chain kinase (MLCK) inhibitor with an IC₅₀ of 50 nM, and inhibits CaM kinase II only at 4000-fold higher concentrations.