
Lucigenin

产品编号: D60110

CAS: 2315-97-1

分子式: C₂₈H₂₂N₄O₆

纯度: 98%

InChi Key: InChIKey=KNJDBYZZKAZQNG-UHFFFAOYSA-N

Smiles: [O-][N+](O.[O-][N+](O-))=O.C[N+](C1=CC=CC=C(C1=CC=CC=C2)C1=CC=CC=C2)C1=CC=CC=C2

外观: 固体粉末

溶解性: DMSO:80mg/mL (156.71 mM)

保存条件: RT

产品介绍: Lucigenin (L-6868) is a chemiluminescent probe used to indicate the presence of endogenously generated superoxide anion radicals in cells.