
Leonurine hydrochloride

产品编号: D60117

CAS: 24735-18-0

分子式: C₁₄H₂₂ClN₃O₅

纯度: 99%

InChi Key: InChIKey=GHVZIMAVDJZQGP-UHFFFAOYSA-N

Smiles: Cl.COC1cc(cc(OC)c1O)C(=O)OCCCCNC(N)=N

外观: 固体粉末

作用通路: Autophagy

作用靶点: Autophagy; TLR

溶解性: DMSO: ≥ 31 mg/mL (89.13 mM) H₂O: 5 mg/mL (14.38 mM; Need ultrasonic)

保存条件: RT

产品介绍: Leonurine hydrochloride is an alkaloid isolated from Leonurus artemisia, with anti-oxidative and anti-inflammatory.