

---

## Docosatrienoic Acid

产品编号: D60126

CAS: 28845-86-5

分子式: C<sub>22</sub>H<sub>38</sub>O<sub>2</sub>

纯度: 98%

外观: 固体粉末

作用靶点: Human Endogenous Metabolite

溶解性: Soluble in Ethanol

保存条件: 2-8°C

产品介绍: Docosatrienoic acid is a rare  $\omega$ -3 fatty acid; inhibits LTB<sub>4</sub> binding to pig neutrophil membranes with an K<sub>i</sub> of 5  $\mu$ M.