
SDMA

产品编号: D60133

CAS: 30344-00-4

分子式: C₈H₁₈N₄O₂

纯度: 98%

InChi Key: InChIKey=HVPFXCBJHIIJGS-LURJTMIESA-N

Smiles: N[C@@H](CCC\N=C(/NC)NC)C(=O)O

外观: 固体粉末

作用通路: Metabolism

作用靶点: Endogenous Metabolite

溶解性: Soluble in DMSO、H₂O

保存条件: 2-8°C

产品介绍: SDMA (Symmetric dimethylarginine) is an endogenous inhibitor of nitric oxide (NO) synthase activity.