
ANGIOTENSIN I

产品编号: D60163

CAS: 484-42-4

分子式: C₆₂H₈₉N₁₇O₁₄

纯度: 98%

InChi Key: InChIKey=ORWYRWWDCYOMK-HBZPZAIKSA-N

Smiles: CC(C)C[C@@H](C(O)=O)NC([C@H](CC1=CNC=N1)NC([C@H](CC2=CC=CC=C2)NC([C@H]3N(C([C@H](CC4=CNC=N4)NC([C@H]([C@@H](C)CC)NC([C@H](CC5=CC=C(O)C=C5)NC([C@H](C(C)C)NC([C@H](CCCNC(N)=N)NC([C@H](CC(O)=O)N)=O)=O)=O)=O)=O)=O)CCC3)=O)=O)=O

外观: 固体粉末

作用靶点: Human Endogenous Metabolite

溶解性: Soluble in H₂O

保存条件: -20°C

产品介绍: Angiotensin I (human, mouse, rat) is the precursor to the vasoconstrictor peptide angiotensin II, cleaved by the angiotensin-converting enzyme (ACE).