
E3 ligase Ligand 2

产品编号: D60165

CAS: 5054-59-1

分子式: C₁₃H₁₀N₂O₅

纯度: 97%

InChi Key: InChIKey=XMPJICVFSDYOEG-UHFFFAOYSA-N

Smiles: Oc1cccc(C(N2C(CCC(N3)=O)C3=O)=O)c1C2=O

外观: 固体粉末

作用靶点: Cereblon; Ligand for E3 Ligase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Thalidomide-4-OH (Cereblon ligand 2) is the Thalidomide-based Cereblon ligand used in the recruitment of CRBN protein. Thalidomide-4-OH (Cereblon ligand 2) can be connected to the ligand for protein by a linker to form PROTACs