

---

## TMB(monosulfate)

产品编号: D60177

CAS: 54827-18-8

分子式: C<sub>16</sub>H<sub>22</sub>N<sub>2</sub>O<sub>4</sub>S

纯度: 99%

外观: 固体粉末

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: TMB monosulfate is a chromogenic substrate used in staining procedures in immunohistochemistry as well as being a visualizing reagent used in enzyme-linked immunosorbent assays (ELISA).