
TG003

产品编号: D60214

CAS: 719277-26-6

分子式: C₁₃H₁₅NO₂S

纯度: 99%

外观: 固体粉末

作用靶点: mCLK4; mCLK1; mCLK2

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: TG003 is a potent Clk1/Sty inhibitor that inhibits Clk1 and Clk4, inhibits cancer cell growth and induces apoptosis, and can be used in the study of Duchenne muscular dystrophy.