
Teneligliptin

产品编号: D60223

CAS: 760937-92-6

分子式: C₂₂H₃₀N₆O₅

纯度: 99%

外观: 固体粉末

作用靶点: DPP4

溶解性: Soluble in DMSO

保存条件: 2-8°C

产品介绍: Teneligliptin (MP-513) is a potent, orally available, competitive, and long-lasting DPP-4 inhibitor. Teneligliptin competitively inhibits human plasma, rat plasma, and human recombinant DPP-4 in vitro, with IC₅₀s of approximately 1 nM