
Cabergoline

产品编号: D60233

CAS: 81409-90-7

分子式: C₂₆H₃₇N₅O₂

纯度: 99%

InChi Key: InChIKey=KORNTPPJEAJQIU-KJXAQDMKSA-N

Smiles: [H][C@@]12Cc3c[nH]c4cccc(c34)[C@@]1([H])C[C@H](CN2CC=C)C(=O)N(CCCN(C)C)C(=O)NC

外观: 固体粉末

作用通路: Autophagy

作用靶点: D3 Receptor; D2 Receptor; Dopamine Receptor; 5-HT Receptor; Autophagy

溶解性: Soluble in DMSO

保存条件: 2-8°C

产品介绍: Cabergoline is an ergot derived-dopamine D2-like receptor agonist that has high affinity for D2, D3, and 5-HT_{2B} receptors (K_i=0.7, 1.5, and 1.2, respectively).