
PSI-6130

产品编号: D60235

CAS: 817204-33-4

分子式: C₁₀H₁₄FN₃O₄

纯度: 99%

InChi Key: InChiKey=NYPiRLYMDJMKGW-VPCXQMTMSA-N

Smiles: C[C@@]1([C@H](N(C=CC(N)=N2)C2=O)O[C@H](CO)[C@H]1O)F

外观: 固体粉末

作用通路: Microbiology/Virology

作用靶点: HCV replication; HCV Protease

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: PSI-6130 is a potent and selective inhibitor of HCV NS5B polymerase, and inhibits HCV replication with a mean IC₅₀ of 0.6 μM.