

---

## 6-Methoxy-N-(3-sulfopropyl)quinolinium

产品编号: D60236

CAS: 83907-40-8

分子式: C<sub>13</sub>H<sub>15</sub>NO<sub>4</sub>S

纯度: 98%

InChi Key: InChIKey=ZVDGOJFPFMINBM-UHFFFAOYSA-N

Smiles: COc1ccc2[n+](CCCS([O-])(=O)=O)cccc2c1

外观: 固体粉末

作用靶点: Fluorescent Dye

溶解性: Soluble in DMSO、H<sub>2</sub>O

保存条件: 2-8°C

产品介绍: SPQ is being used to examine and measure membrane chloride transport mechanisms.