
FMOC-3VVD-OH

产品编号: D60247

CAS: 863971-44-2

分子式: C₃₆H₅₁N₃O₇

纯度: 98%

InChi Key: InChIKey=SMVBIYSFDNRJJI-FCUBIXSHSA-N

Smiles: C(OC(N([C@H](C(N[C@H](C(N([C@H]([C@@H](CC(O)=O)OC)[C@H](CC)C)C)=O)[C@@H](C)C)=O)C(C)C)=O)C1C=2C(C=3C1=CC=CC3)=CC=CC2

外观: 固体粉末

作用靶点: Cleavable

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Fmoc-3VVD-OH is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs)