

---

## Thalidomide-O-COOH

产品编号: D60279

CAS: 1061605-21-7

分子式: C<sub>15</sub>H<sub>12</sub>N<sub>2</sub>O<sub>7</sub>

纯度: 96%

InChi Key: InChIKey=ZVWIXIGBWIEDFO-UHFFFAOYSA-N

Smiles: O=C(N(C(CCC1=O)C(N1)=O)C2=O)C(C2=CC=C3)=C3OCC(O)=O

外观: 固体粉末

作用靶点: Cereblon; Ligand for E3 Ligase

溶解性: Soluble in DMSO

保存条件: 2-8°C

产品介绍: Thalidomide-O-COOH (Cereblon ligand 3) is the Thalidomide-based Cereblon ligand used in the recruitment of CRBN protein. Thalidomide-O-COOH (Cereblon ligand 3) can be connected to the ligand for protein by a linker to form PROTACs