
MW-150

产品编号: D60297

CAS: 1628502-91-9

分子式: C₂₄H₂₃N₅

纯度: 98%

InChi Key: InChIKey=CIIVUDIZZJLXCN-UHFFFAOYSA-N

Smiles: CN1CCN(CC1)c1cc(-c2ccncc2)c(nn1)-c1ccc2ccccc2c1

外观: 固体粉末

作用通路: Autophagy

作用靶点: p38 α ; p38MAPK; Autophagy

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: MW150 (MW01-18-150SRM) is a selective, CNS penetrant, and orally active inhibitor of p38 α MAPK with a K_i of 101 nM. MW-150 inhibits the ability of the endogenous p38 α MAPK to phosphorylate an endogenous substrate MK2 in activated glia