

---

## Boc-Val-Cit-PAB-PNP

产品编号: D60305

CAS: 870487-10-8

分子式: C<sub>30</sub>H<sub>40</sub>N<sub>6</sub>O<sub>10</sub>

纯度: 98%

InChi Key: InChIKey=SKWLDDHVHXQPOG-ZEQLZLVSA-N

Smiles: CC(C)[C@H](NC(=O)OC(C)(C)C(=O)N[C@@H](CCNC(N)=O)C(=O)Nc1ccc(COC(=O)Oc2ccc(cc2)[N+](=[O-])=O)cc1

外观: 固体粉末

作用靶点: Protease Cleavable; Cleavable

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Boc-Val-Cit-PAB-PNP is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs)