
AMG 208

产品编号: D50653

CAS: 1002304-34-8

保存条件: Store at -20°C for 3 year(Powder); In DMSO or others solvent store at 2-4°C for two weeks, at -20°C for six months.

产品介绍: CAS: 11002304-34-8

分子式: C₂₂H₁₇N₅O₂

分子量: 383. 41

MDL: MFCD17215205

纯度: ≥99% (HPLC)

溶解性: 溶于DMSO

作用靶点: c-Met/HGFR;

作用通路: Protein Tyrosine Kinase/RTK;

产品描述: AMG-208是c-Met抑制剂, IC₅₀为9.3 nM。