
BI-4464

产品编号: D60362

CAS: 1227948-02-8

分子式: C₂₈H₂₈F₃N₅O₄

纯度: 98%

InChi Key: InChIKey=QUSSZSMDFABHLI-UHFFFAOYSA-N

Smiles: COc1cc(ccc1Nc1ncc(c(Oc2cccc3CCC(=O)c23)n1)C(F)(F)F)C(=O)NC1CCN(C)CC1

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: PTK2/FAK; FAK; Ligands for Target Protein for PROTAC

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: BI-4464 is a highly selective ATP competitive inhibitor of PTK2/FAK, with an IC₅₀ of 17 nM. A PTK2 ligand for PROTAC