
AMG-333

产品编号: D60368

CAS: 1416799-28-4

分子式: C₂₀H₁₂F₅N₃O₄

纯度: 98%

InChi Key: InChIKey=QEBYISWYMFIXOZ-INIZCTEOSA-N

Smiles: O=C(O)c1ccc(nc1)C(=O)N[C@H](c2ncccc2F)c3ccc(OC(F)(F)F)c(F)c3

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: TRPM8; TRP/TRPV Channel

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: AMG 333 is TRPM8 antagonist with an IC₅₀ of 13 nM, a Clinical Candidate for the Treatment of Migraine.