
UCB9608

产品编号: D60374

CAS: 1616413-96-7

分子式: C₂₀H₂₆N₈O₂

纯度: 98%

InChi Key: InChIKey=WRONAJQPZWDYAR-ZDUSSCGKSA-N

Smiles: C[C@H]1CN(CCN1C2=NC(=NC3=C2C=NN3C)N)C(=O)NC4=C(C=C(C=C4)OC)C

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PI4K

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: UCB9608 is a potent, selective and orally active PI4KIII β inhibitor, with an IC₅₀ of 11 nM, selective over PI3KC2 α , β , and γ lipid kinases. UCB9608 improves metabolic stability and exhibits excellent pharmacokinetic profile, acts as a potent immunosuppressive agent