
SLV-2436

产品编号: D60387

CAS: 2095704-43-9

分子式: C₁₉H₁₅ClN₄O

纯度: 98%

InChi Key: InChIKey=YQVUADHJKWJHAF-UHFFFAOYSA-N

Smiles: O=C1C=CC(=CN1CC=2C=CC=C(Cl)C2)C=3C=CC=4C(=NNC4C3)N

外观: 固体粉末

作用通路: MAPK

作用靶点: MNK2; MNK1; MNK

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: SLV-2436 is a highly potent and ATP-competitive inhibitor of MNK1 and MNK2 with IC₅₀s of 10.8 nM and 5.4 nM, respectively.