
NH2-PEG5-OH

产品编号: D60394

CAS: 34188-11-9

分子式: C₁₀H₂₃NO₅

纯度: 97%

外观: 固体粉末

作用靶点: PEGs; Non-cleavable

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: NH₂-PEG5-OH is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs.

NH₂-PEG5-OH is also a non-cleavable 5 unit PEG ADC linker used in the synthesis of antibody-drug conjugates (ADCs)