
Homo-PROTAC pVHL30 degrader 1

产品编号: D60427

CAS: 2244684-49-7

分子式: C₅₈H₈₂N₈O₁₄S₂

纯度: 98%

InChi Key: InChIKey=WGJCHHJGGFCCRS-DEYDLUNASA-N

Smiles: O=C(N[C@@H](C(C)(C)C)C(N1[C@H](C(NCC2=CC=C(C3=C(C)N=CS3)C=C2)=O)C[C@@H](O)C1)=O)COCCOCCOCCOCCOCCOCC(N[C@@H](C(C)(C)C)C(N4[C@H](C(NCC5=CC=C(C6=C(C)N=CS6)C=C5)=O)C[C@@H](O)C4)=O)=O

外观: 固体粉末

作用靶点: Ligand for E3 Ligase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Homo-PROTAC pVHL30 degrader 1 is a potent pVHL30 degrader based on PROTAC.