
Uprosertib

产品编号: D60568

CAS: 1047634-65-0

分子式: C₁₈H₁₆Cl₂F₂N₄O₂

纯度: 99%

InChi Key: InChIKey=AXTAPYRUEKNRBA-JTQLQIEISA-N

Smiles: Cn1c(c(c1)Cl)c1c(oc(c1)C(=O)N[C@@H](Cc1cc(c(cc1)F)F)CN)Cl

外观: 固体粉末

作用通路: Cytoskeletal Signaling

作用靶点: Akt

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: Uprosertib (GSK2141795) is a potent, selective Akt broad-spectrum inhibitor with IC₅₀ values of 180/328/38 nM for Akt1/Akt2/Akt3.