
LY-2584702 tosylate salt

产品编号: D60580

CAS: 1082949-68-5

分子式: C₂₈H₂₇F₄N₇O₃S

纯度: 98%

InChi Key: InChIKey=HDYUXDNMHBQKAU-UHFFFAOYSA-N

Smiles: Cc1ccc(cc1)S(=O)(=O)O.Cn1cc(nc1C1CCN(CC1)c1ncnc2c1cn[nH]2)c1cc(c(cc1)F)C(F)(F)F

外观: 固体粉末

作用通路: MAPK

作用靶点: S6 Kinase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: LY-2584702 tosylate salt is a selective ATP competitive inhibitor of p70S6K with an IC₅₀ of 4 nM. In S6K1 enzyme assay, the IC₅₀ of LY-2584702 is 2 nM