
AZ3146

产品编号: D60590

CAS: 1124329-14-1

分子式: C₂₄H₃₂N₆O₃

纯度: 99%

InChi Key: InChIKey=YUKWVHPTFRQHMF-UHFFFAOYSA-N

Smiles: COc1c(Nc2ncc3n(C)c(=O)n(C4CCCC4)c3n2)ccc(OC2CCN(C)CC2)c1

外观: 固体粉末

作用通路: Cytoskeletal Signaling

作用靶点: Kinesin

溶解性: DMSO: 100 mg/mL (220.97 mM; Need ultrasonic)

保存条件: -20°C

产品介绍: AZ3146 is a reasonably potent Mps1 and TTK inhibitor, with IC₅₀ of 35 nM for Mps1Cat