
SNS-314 Mesylate

产品编号: D60595

CAS: 1146618-41-8

分子式: C₁₈H₁₅ClN₆O₅S₂·CH₄O₃S

纯度: 99%

InChi Key: InChIKey=FYCOPVDEFFWSR-UHFFFAOYSA-N

Smiles: CS(=O)(=O)O.Clc1cccc(NC(=O)Nc2ncc(CCNc3ncnc4c3scc4)s2)c1

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint

作用靶点: Aurora Kinase

溶解性: DMSO: 150 mg/mL (284.61 mM; Need ultrasonic)

保存条件: -20°C

产品介绍: SNS-314 mesylate is a potent and selective aurora kinase inhibitor with IC₅₀s of 9, 31, and 6 nM for aurora A, B and C, respectively