
Sal003

产品编号: D60599

CAS: 1164470-53-4

分子式: C₁₉H₁₆Cl₄N₂O₂S

纯度: 98%

InChi Key: InChIKey=TVNBASWNLOIQML-IZZDOVSWSA-N

Smiles: ClC1=CC=C(NC(=S)NC(NC(=O)\C=C\C2=CC=CC=C2)C(Cl)(Cl)Cl)C=C1

外观: 固体粉末

作用靶点: Apoptosis; PERK; Phosphatase

溶解性: DMSO: ≥ 100 mg/mL (215.88 mM)

保存条件: 2-8°C

产品介绍: Sal003 is a potent, specific and cell-permeable inhibitor of the eukaryotic translation initiation factor 2α (eIF2α) phosphatase. Sal003 is a derivative of salubrinal