
AP-III-a4

产品编号: D60601

CAS: 1177827-73-4

分子式: C₃₁H₄₃FN₈O₃

纯度: 96%

InChi Key: InChiKey=MOVYITHKOHMLHC-UHFFFAOYSA-N

Smiles: C1CCC(CC1)CNc1nc(nc(n1)NCc1ccc(cc1)F)Nc1ccc(cc1)CC(=O)NCCOCCOCCN

外观: 固体粉末

作用通路: Metabolism

作用靶点: Glucokinase

溶解性: Soluble in DMSO、Ethanol、H₂O

保存条件: -20°C

产品介绍: AP-III-a4 (ENOblock) is the first nonsubstrate analogue inhibitor of enolase with IC₅₀ of 0.576 μM.