
GSK2334470

产品编号: D60618

CAS: 1227911-45-6

分子式: C₂₅H₃₄N₈O

纯度: 99%

InChi Key: InChIKey=QLPHOXTXAKOFMU-WBVHZDCISA-N

Smiles: CNc1nc(cc(n1)-c1ccc2c(N)n[nH]c2c1)N1C[C@H](CC[C@H]1C)C(=O)NC1CCCCC1

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PDK

溶解性: DMSO: ≥ 50 mg/mL (108.09 mM)

保存条件: -20°C

产品介绍: GSK2334470 is a highly specific and potent inhibitor of PDK1 with an IC₅₀ of 10 nM.