
XMD8-87

产品编号: D60628

CAS: 1234480-46-6

分子式: C₂₄H₂₇N₇O₂

纯度: 98%

InChi Key: InChiKey=LGLHCXISMKHLIK-UHFFFAOYSA-N

Smiles: COc1cc(ccc1Nc1ncc2NC(=O)c3ccccc3N(C)c2n1)N1CCN(C)CC1

外观: 固体粉末

作用通路: DNA Damage/DNA Repair

作用靶点: PPAR; Tyrosinase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: XMD8-87 is a potent TNK2 inhibitor with IC₅₀ values of 38 and 113 nM for the D163E and R806Q mutations, respectively