
XMD16-5

产品编号: D60659

CAS: 1345098-78-3

分子式: C₂₃H₂₄N₆O₂

纯度: 98%

InChi Key: InChiKey=AGLKBEPKKDHHKY-UHFFFAOYSA-N

Smiles: CN1c2cccc2C(=O)Nc2cnc(Nc3ccc(cc3)N3CCC(O)CC3)nc12

外观: 固体粉末

作用通路: DNA Damage/DNA Repair

作用靶点: PPAR; Tyrosinase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: XMD16-5 is a potent TNK2 inhibitor with IC₅₀ values of 16 and 77 nM for the D163E and R806Q mutations, respectively.