
G007-LK

产品编号: D60668

CAS: 1380672-07-0

分子式: C₂₅H₁₆ClN₇O₃S

纯度: 99%

InChi Key: InChIKey=HIWVLHPKZNBSBE-OUKQBFOZSA-N

Smiles: CS(=O)(=O)c1cnc(cc1)c1nnc(n1c1ccccc1Cl)/C=C/c1nnc(o1)c1ccc(cc1)C#N

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: PARP; Wnt/beta-catenin

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: G007-LK is a selective inhibitor of TNKS1 and TNKS2, with IC₅₀s of 46 nM and 25 nM, respectively.