

---

## ZM241385

产品编号: D60671

CAS: 139180-30-6

分子式: C<sub>16</sub>H<sub>15</sub>N<sub>7</sub>O<sub>2</sub>

纯度: 96%

InChi Key: InChIKey=PWTBZOIUWZOPFT-UHFFFAOYSA-N

Smiles: c1(ccc(cc1)CCNc1nc2n(c(n1)N)nc(n2)c1occc1)O

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: Adenosine Receptor

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: ZM-241385 is a high-affinity antagonist ligand selective for the adenosine A<sub>2A</sub> receptor.