
XEN445

产品编号: D60691

CAS: 1515856-92-4

分子式: C₁₈H₁₇F₃N₂O₃

纯度: 96%

InChi Key: InChIKey=NBGRERFNOKZQLO-AWEZLNQCLSA-N

Smiles: c1(ccc(c(c1)C(=O)O)N1C[C@H](CC1)OCc1ncccc1)C(F)(F)F

外观: 固体粉末

作用通路: Metabolism

作用靶点: Lipase

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: XEN445, a potent and selective endothelial lipase(EL) inhibitor (IC₅₀=0.237 uM), exhibitis good ADME and PK properties.