

---

## SCH 58261

产品编号: D60700

CAS: 160098-96-4

分子式: C<sub>18</sub>H<sub>15</sub>N<sub>7</sub>O

纯度: 99%

InChi Key: InChIKey=UTLPKQYUXOEJIL-UHFFFAOYSA-N

Smiles: Nc1nc2n(CCc3ccccc3)ncc2c2nc(nn12)-c1ccco1

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: Adenosine Receptor

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: SCH 58261 is a potent, selective and competitive antagonist of adenosine A<sub>2A</sub> receptor with an IC<sub>50</sub> of 15 nM, and displays 323-, 53- and 100-fold more selective for A<sub>2A</sub> receptor than A<sub>1</sub>, A<sub>2B</sub>, and A<sub>3</sub> receptors, respectively