
LJI308

产品编号: D60704

CAS: 1627709-94-7

分子式: C₂₁H₁₈F₂N₂O₂

纯度: 99%

InChi Key: InChIKey=YUYJEQHNWKQNB-TUHFFFAOYSA-N

Smiles: Oc1c(F)cc(cc1F)-c1cnccc1-c1ccc(cc1)N1CCOCC1

外观: 固体粉末

作用通路: MAPK

作用靶点: S6 Kinase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: LJI308 is a potent, and pan-RSK (p90 ribosomal S6 kinase) inhibitor with IC₅₀ of 6 nM, 4 nM, and 13 nM for RSK1, RSK2, and RSK3, respectively.