
EAI045

产品编号: D60735

CAS: 1942114-09-1

分子式: C₁₉H₁₄FN₃O₃S

纯度: 99%

InChi Key: InChIKey=YTUFHOKUFOQRDF-UHFFFAOYSA-N

Smiles: c1ccc2c(c1)CN(C2=O)C(C(=O)Nc1scn1)c1cc(ccc1O)F

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: EGFR

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: EAI045, an allosteric inhibitor, targets towards drug-resistant EGFR mutants but avoids the wild-type receptor.