
SL327

产品编号: D60780

CAS: 305350-87-2

分子式: C₁₆H₁₂F₃N₃S

纯度: 98%

InChi Key: InChIKey=JLOXTZFYJNCPIS-UHFFFAOYSA-N

Smiles: NC(Sc1ccc(N)cc1)=C(C#N)c1ccccc1C(F)(F)F

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint

作用靶点: DNA/RNA Synthesis; MEK

溶解性: DMSO: 68 mg/mL (202.77 mM; Need ultrasonic) Ethanol: 0.1 mg/mL (0.30 mM; Need ultrasonic and warming)

保存条件: -20°C

产品介绍: SL327 inhibits MEK1 and MEK2, with IC₅₀ values of 180 nM and 220 nM, respectively.