
T0070907

产品编号: D60784

CAS: 313516-66-4

分子式: C₁₂H₈ClN₃O₃

纯度: 99%

InChi Key: InChIKey=FRPJSHKMZHWJBE-UHFFFAOYSA-N

Smiles: c1cc(c(cc1[N+](=O)[O-])C(=O)Nc1ccncc1)Cl

外观: 固体粉末

作用通路: DNA Damage/DNA Repair

作用靶点: PPAR

溶解性: DMSO: 62.5 mg/mL (225.10 mM; ultrasonic and warming and heat to 60°C) H₂O: 1.1 mg/mL (3.96 mM; Need ultrasonic)

保存条件: -20°C

产品介绍: T0070907 is a potent PPAR γ antagonist with a K_i of 1 nM.