
SKF 38393 hydrochloride

产品编号: D60856

CAS: 62717-42-4

分子式: C₁₆H₁₇NO₂

纯度: 99%

InChi Key: InChIKey=YEWBJCLOUYPAOH-UHFFFAOYSA-N

Smiles: Cl.OC1=C(O)C=C2C(CNCCC2=C1)C1=CC=CC=C1

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: Dopamine Receptor

溶解性: DMSO: ≥ 34 mg/mL (116.53 mM) H₂O: 10 mg/mL (34.27 mM; Need ultrasonic)

保存条件: RT

产品介绍: SKF 38393 hydrochloride is a selective agonist of the dopamine D1 receptor (D1DR) with an IC₅₀ of 110 Nm