

---

## Givinostat hydrochloride monohydrate

产品编号: D60882

CAS: 732302-99-7

分子式: C<sub>24</sub>H<sub>30</sub>ClN<sub>3</sub>O<sub>5</sub>

纯度: 98%

InChi Key: InChIKey=FKGKZBBDJSKCIS-UHFFFAOYSA-N

Smiles: CCN(CC)CC1=CC2=C(C=C1)C=C(C=C2)COC(=O)NC3=CC=C(C=C3)C(=O)NO.O.Cl

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: HDAC

溶解性: Soluble in DMSO、H<sub>2</sub>O

保存条件: -20°C

产品介绍: Givinostat hydrochloride monohydrate (ITF-2357 hydrochloride monohydrate) is a HDAC inhibitor with an IC<sub>50</sub> of 198 and 157 nM for HDAC1 and HDAC3, respectively.